



GAMMAVITA® Phase/Suspension Mixer

For processing liquid to semi-solid products and as a melting tank

GAMMAVITA® Phase/Suspension Mixer Quality in every detail

For economic production in a highly competitive environment, you need the highest level of quality and reliability, as provided by the NETZSCH Vakumix *GAMMA-VITA*® series phase tanks.

The design of these machines guarantees repeatable production results with maximum efficiency for processing liquid to semi-solid products and as a melting tank.

- Melting
- Mixing
- Homogenization
- Dispersing
- Degassing



Melting of wax pellets



Melting of wax blocks

Patented Circulating Agitator with grinding head





Mixing inside the tank



Mixing via circulation line and induction of powder

GAMMAVITA® Phase/Suspension Mixer Technical Features

- Pressure- and Vacuum-tight design
- GMP-compliant machine design
- Scalability
- CIP- and, optionally, SIP-capable machine design
- All product-wetted parts in high-quality stainless steel (316L)
- Short melting and mixing times due to external circulation line
- Patented (EP 0760254 A1) Circulating Agitator:
 - Optional recirculation of product inside the tank or via external circulation line
 - Installation at the lowest point of the tank bottom
 - Stepless control of turbulence intensity via frequency-controlled drive
 - No additional pumps necessary for CIP and product transfer

- Virtually complete emptying of the machine via the Circulating Agitator
- Additional measuring probe assembly in the housing of the Circulation Agitator possible (temperature, pH value, pressure, etc.)
- Can be optionally equipped with homogenizer instead of Circulation Agitator
- Various suction options for solids and liquids in a closed system
- Can be heated/cooled via energy modules with features and media to meet customer requirements
- Vacuum operation via water ring vacuum pumps
- Various machine control and automation styles available
- Can be combined with the model series KAPPAVITA®, DELTAVITA® and OMEGA® VITA as SYSTEMVITA
- Sizes from 10 12,000 liter (larger sizes upon inquiry)



System VITA®, consisting of a DELTAVITA® 2000 agitator bead mill and a GAMMAVITA® 60 batch tank



KAPPAVITA® HM 500 homogenizing mixer with GAMMAVITA® 500 phase tank

GAMMAVITA® Phase/Suspension Mixer Models and Sizes

Available sizes and executions:

- Sizes from 10 12.000 liter
- *GAMMAVITA*[®] with circulation agitator
- GAMMAVITA® with rotor/stator homogenizer
- GAMMAVITA® with magnetical stirrer

All GAMMA VITA® models are optionally available with top driven stirrers, like the NETZSCH Vakumix Easy Visco.



GAMMAVITA® 275 – mobile with work platform



GAMMAVITA® 140 – with Load cells

GAMMAVITA® 300 – mobile

The NETZSCH Group is an owner-managed, international technology company with headquarters in Germany. The Business Units Analyzing & Testing, Grinding & Dispersing and Pumps & Systems represent customized solutions at the highest level. More than 3,700 employees in 36 countries and a worldwide sales and service network ensure customer proximity and competent service.

Our performance standards are high. We promise our customers Proven Excellence – exceptional performance in everything we do, proven time and again since 1873.

Proven Excellence.

Business Unit Grinding & Dispersing – The World's Leading Grinding Technology

NETZSCH-Feinmahltechnik – Germany NETZSCH Trockenmahltechnik – Germany NETZSCH Vakumix – Germany NETZSCH Lohnmahltechnik – Germany NETZSCH Mastermix – Great Britain NETZSCH FRÈRES – France NETZSCH España – Spain ECUTEC – Spain NETZSCH Machinery and Instruments – China NETZSCH Technologies India Private – India NETZSCH Tula – Russia NETZSCH Makine Sanayi ve Ticaret – Turkey NETZSCH Korea – Korea NETZSCH Premier Technologies – USA NETZSCH Equipamentos de Moagem – Brazil

NETZSCH Vakumix GmbH Zeppelinstraße 1 28844 Weyhe-Dreye Germany Tel.: +49 421 8 49 56 0 Fax: +49 421 8 49 56 90 info.nvm@netzsch.com

