




OmegaVita® High-pressure Homogenizer

For the production of emulsions with droplet sizes in the nanometer range

OmegaVita® High-pressure Homogenizer Dispersion of liquid to semi-solid suspensions


Quality in every detail

For economic production in a highly competitive environment, you need the highest level of quality and reliability, as provided by the NETZSCH Vakumix *OMEGA-VITA®* series high-pressure homogenizers. The design of these machines guarantees repeatable production results with maximum efficiency. A typical application is the production of nanoemulsions.



Technical Features

- High-pressure homogenizer up to 700 bar
- Efficient dispersion with the Netzsch Dispersion -Device (NDD) a combination of nozzle and valve technology
- Batch sizes from 15 milliliters to 10000 liters.
- Scalable
- SPS control
- Operator management with password protection in different levels
- Optional formulation management and data logging in accordance with CFR part 11/21
- GMP-compliant machine design
- CIP- and, optionally, SIP-capable machine design
- All product-wetted parts are manufactured and certified according to FDA standards
- Material and calibration certificates are delivered with the machine
- Comprehensive documentation in the national language, FAT, IQ and OQ documents
- Can be combined with the model series *KAPPAVITA*® and *GAMMAVITA*® as System *VITA*



and emulsions

OmegaVita® series sizes

		ОмедаLав	ОмеgaVita® 60	ОмеgaVita® 500	OMEGAVITA® 2000
Pressure	[bar]	10 - 700	10 - 700	10 - 700	10 - 700
Batch sizes (min max.)	[1]	0,3 - 3	5 - 100	20 - 2 000	100 - 10 000
Throughput	[1]	1 - 20	18 - 60	20 - 500	600 - 2000



The NETZSCH Group is an owner-managed, international technology company with headquarters in Germany. The Business Units Analyzing & Testing, Grinding & Dispersing and Pumps & Systems represent customized solutions at the highest level. More than 3,700 employees in 36 countries and a worldwide sales and service network ensure customer proximity and competent service.

Our performance standards are high. We promise our customers Proven Excellence – exceptional performance in everything we do, proven time and again since 1873.

Proven Excellence.

Business Unit Grinding & Dispersing – The World's Leading Grinding Technology

NETZSCH-Feinmahltechnik – Germany

NETZSCH Trockenmahltechnik – Germany

NETZSCH Vakumix – Germany

NETZSCH Lohnmahltechnik – Germany

NETZSCH Mastermix – Great Britain

NETZSCH FRÈRES – France

NETZSCH España – Spain

ECUTEC - Spain

NETZSCH Machinery and Instruments – China

NETZSCH Technologies India Private – India NETZSCH Tula – Russia

NETZSCH Makine Sanayi ve Ticaret – Turkey

NETZSCH Korea – Korea

NETZSCH Premier Technologies – USA

NETZSCH Equipamentos de Moagem – Brazil

NETZSCH Vakumix GmbH Zeppelinstraße 1 28844 Weyhe-Dreye Germany

Tel.: +49 421 8 49 56 0 Fax: +49 421 8 49 56 90 info.nvm@netzsch.com

